



cantabria labs

HELIOCARE

LUMINANCE

THE ULTIMATE SOLUTION FOR BRIGHT, EVEN SKIN TONE



Fernblock
Researched for
more than 12 years in




CLINICAL EVIDENCE ON HELIOCARE LUMINANCE

In-vivo clinical study on the effect of Heliocare Luminance on hyperpigmentation in Asian women.

Method of study:

- This clinical study is randomized and placebo-controlled.
- 26 subjects and 12 subjects (Females with spots, uneven complexion tone and aging signs) were assigned to the Active Luminance group and the Placebo group, respectively.
- Evaluation: Baseline, Week 4, Week 8, and Week 12.


The purpose of this study was to assess the efficacy of Heliocare Luminance to prevent hyperpigmentation and skin damages induced by solar radiation and to improve area, length and depth of wrinkles, skin firmness and elasticity, skin barrier function and skin hydration in female subjects.






Global improvement of Area, Length and Depth of wrinkles

After 12 weeks of study, it was demonstrated that Heliocare Luminance has efficacy to improve hyperpigmentation, skin elasticity and firmness, deep wrinkles, skin barrier function and skin hydration.

WEEK 0



WEEK 12



Data on file, Eurofins Korea, 2023.



HELIOCARE LUMINANCE

is recommended for the following indications:

- Uneven skin complexion and hyperpigmentation
- Post dermatological procedures to help prevent formation of spots and hyperpigmentation
- Adjuvant to aesthetic treatments or lightening procedure, to receive optimal results
- Dehydrated skin with impaired skin barrier function
- Diminish appearance of fine lines and wrinkles, achieve a youthful looking skin

Dosage

Take 2 capsules daily, preferably one in the morning and one at night.

Contains 60 capsules, sufficient to cover the natural 28-day skin renewal cycle.

Active ingredients of each capsule

- 240mg Fernblock®
- 125mg L-Glutathione
- 20mg Rosemary Extract
- 8mg Nicotinamide (Vitamin B3)



Brought to you by:

NEOASIA

NEOASIA (S) PTE LTD Co. Reg. No. 199504326C
6 Tagore Drive #04-07 Tagore Building Singapore 787623
+65 6552 7787 | info@neoasia.com.sg | www.neoasia.com.sg


cantabria labs
celebrate life

SUPPLEMENTS DIFFERENTIATION CHART



	Heliocare Fern Mar Plus/ D Plus	Heliocare C Plus Oral Capsules	Heliocare Fern PLE+	Heliocare 60 Oral Capsules	Heliocare Luminance
Where to buy	Clinics exclusive	Clinics exclusive	Clinics exclusive	Retail & Online stores	Clinics exclusive
Key ingredient	Fernblock® 480mg	Fernblock®+ 240mg	Fernblock® 240mg	Fernblock® 240mg	Fernblock® 240mg
Additional ingredients	Lycopene Vitamin C & E Vitamin D	Prebiotics Vitamins B3, C & E Cystine	Pomegranate Extract Vitamin C	-	L-glutathione Rosemary Extract Niacinamide
Who should take	Intense sun exposure, sports, sun sensitive skin and post dermatological care.	Those that are mostly indoor and exposed to blue light such as mobile devices, laptop and television. Additional protection against infrared radiation A such as heat from the sun, thermal heat and cooking.	Intense efficacy to target dark spots, uneven skin tone and provide anti-aging benefits. Improve overall skin radiance and glow.	Enhance daily sun protection, beauty, travel, sports.	The ultimate solution for bright, even skin tone. Reduce hyperpigmentation, enhance skin radiance and improve overall skin health.
Capsules per bottle	30 caps	30 caps	60 caps	60 caps	60 caps
Direction of use	Take 1 capsule 30 minutes before sun exposure. If prolonged exposure, retake 3 hours later.	Indoor: Take 1 capsule daily. Outdoor: Take 1 capsule 30 minutes before sun exposure. If prolonged exposure, retake 3 hours later.	Take 2 capsules daily, one in the morning and one at night	For regular sun exposure, take 1 capsule 30 minutes before sun exposure. If prolonged exposure, retake 3 hours later. For intense exposure, take 2 capsules instead.	Take 2 capsules daily, one in the morning and one at night



HELIOCARE

HELIOCARE

360°

Oil-Free Compact

Tinted photoimmunoprotection for
an even skin tone



Broad spectrum
(SPF50+ ; PA ++++; HEV-VIS & IR-A)



Prevents and repair sun damage



Helps prevent hyperpigmentation



Helps conceal blemishes



3 shades: pearl, beige and bronze

INDICATIONS

- Tinted photoprotection that evens skin tone
- Suitable for all skin type including sensitive





HELIOCARE

360°

SPORT

Transparent Stick

High facial photoimmunoprotection
for outdoor sports



Broad spectrum photoimmunoprotection
(SPF50+, PA +++++, Visible and IR)



Very easy to apply, doesn't need to be
spread with hands



Water and sweat resistant.

INDICATIONS

- Sun protection specially formulated for outdoor sports
- On-the-go format





HELIOCARE

360°

Sensation

Universal daily photoimmunoprotection
with a great user experience



Broad spectrum photoimmunoprotection
(SPF50+, PA +++++, Visible & IR-A)



Prevents and repairs sun damage



Ultra light oil-free texture



Great sensorial experience

INDICATIONS



- Pleasant, easy-to-use daily facial photoprotection
- Suitable for all skin types

HELIOCARE

360°

Water Gel

360° photoimmunoprotection and long-lasting hydration

360° DERMATOLOGICAL PHOTOIMMUNOPROTECTION



Full spectrum

(UVA • UVB SPF 50+ • HEVL • IR-A)



Innovative continuous-release hydration system



Ultra-light water gel texture



Eye tolerance and fast absorption

INDICATIONS



- Dry, dehydrated skin
- Sun-related pathologies, photosensitivity, photoaging and spots
- Suitable for all skin types including combination and oily

HELIOCARE

360°

Mineral Tolerance Fluid

360° photoimmunoprotection 100% mineral
with fluid and invisible texture



Full spectrum
(UVA • UVB SPF 50 • HEVL • IR)



100% mineral filter formulation



No whitening effect



Excellent skin tolerance

INDICATIONS



- Intolerant and sensitive skin
- Sun related pathologies, photosensitivity, photoaging and hyperpigmentation
- Protection of the sensitized skin after aesthetic procedures

HELIOCARE

360°

Pigment Solution Fluid

Prevention and correction of hyperpigmentation and dark spots caused by solar radiation



Very high broad-spectrum photoimmunoprotection (UVA • UVB SPF50+ • HEVL • IRA)



Helps to prevent and correct hyperpigmentation and oxidative damage



Light and silky texture that unifies skin tone

INDICATIONS



- Specially indicated for hyperpigmented or prone to spots skin
- Excellent as make-up base

HELIOCARE

360°

Gel Oil-Free

360° photoimmunoprotection with a light dry touch skin feel



Full spectrum
(UVA • UVB SPF 50 • HEVL • IR)



Light oil-free gel texture



Matte effect and dry touch



Sebum-regulating activity

INDICATIONS



- Specifically formulated for combination, oily or acne-prone skin
- Hot and humid climates
- Sun related pathologies, photosensitivity, photoaging and hyperpigmentation